



### Dr. Adjustment Kit

1 green mounted stone, 1 lg brown rubber point, 1 sm brown rubber point, 1 carbide bur, 1 carving knife, sample piece of blue hi-shine paste, 5 mini rag buffs, 5 mini fuzz buffs.  
**4014-02**



### Sample Models

Unbreakable™ partial      **4011-01 ea.**  
iFlex™ partial              **4009-01 ea.**  
Karadent™ full denture    **4016-01 ea.**



### Green Mounted Stones

Ideal for chairside adjustments.  
**4501-01 box of 100**



### Carbide Bur

Specifically designed to cut through thermoplastic materials. These burs are long lasting with high cutting efficiency.  
**50C79-01 ea.**



### Brown Rubber Points

Ideal for chairside adjustments. Use after green mounted stone to smoothen area when adjusting TCS appliances.  
Small **4511-01S ea.** | Large **4510-01L ea.**



### Carving Knife & Replacement Blades

Ideal for removing fibers that sometimes appear on thermoplastic restorations during adjustments. Replacement blades are specific to this knife.  
Knife **4201-01 ea.** | Blades **4202-01 pack of 12**



### Mini Fuzz Buff

Recommended to smoothen TCS appliances after adjustments.  
**4222-01 ea.**



### Mini Rag Buff

Recommended for touch-up polishing of TCS appliances after using the mini fuzz buff. Use with blue hi-shine.  
**4221-01 ea.**



### Blue Hi-Shine

Recommended for touch-up polishing of TCS appliances. Use with the mini rag buff.  
**4215-01 (6.7 oz) ea.**



### ORDER OPTIONS

866-426-2970 | [www.tcsdentalinc.com](http://www.tcsdentalinc.com)