



TrueDent®

One Piece. One Print One Better Denture.



Scan



Send



Deliver

True Aesthetics. Efficient Practice Flow Made Possible.

- **Efficient** – Your try-in is a finished denture.
- **Monolithic** — One seamless print — no separate base and teeth to assemble.
- **Lifelike Color** — Five resin colors blended for natural gradient shading.
- **Precision** — PolyJet accuracy other digital dentures can't match.

TrueDent® Dentures

The technology behind every Peterson denture.

Most digital dentures are built in two pieces — a separate base and teeth, bonded together — with flat, uniform shading. TrueDent® prints the whole denture as one seamless piece, blending five resin colors for natural gradient shading no assembled denture can match.

Turnaround Times:

Standard: 5 days in lab. Expedited options are available on select TrueDent® products.

Frequently Asked Questions

What products does TrueDent® cover?

Full dentures (Immediate, Complete, and Copy Dentures) and Partial dentures (both full acrylic and metal-based).

Can a TrueDent Denture be relined?

Yes. TrueDent™ dentures can be relined using conventional denture reline procedures when clinically indicated. Because TrueDent™ is a printed resin material, an appropriate bonding agent/adhesion protocol should be used to achieve optimal bond strength between the denture base and the reline material.

How is this different from Peterson's TruColor dentures?

It's not a different product — it's the same denture Peterson has always made, now proudly carrying its official Stratasys name, TrueDent®.

Does the lab store a digital copy?

Yes — same as Peterson's existing digital denture workflow, a digital copy is retained.

READY TO LEARN MORE?

Contact your Peterson account manager or sales manager to get started!