



## Need a copy?

Sometimes, making a new smile doesn't mean starting from scratch. Simply scan the patient's existing denture, and we'll create a new one that's stronger than the original.

### You might need a copy if:

- The patient likes their existing denture, but complains about fit.
- The patient just wants to replace what they already have.
- The patient wants a backup.
- You have an upcoming conversion, but the old denture is brittle.
- If you're sending a repair to the lab, we can fix it AND make a back up!


**Let the digital workflow work for you, and get your patients smiling again with a Copy Denture!**



Scan



Send



Deliver

# Copy Dentures

When all you need is an exact copy

## Turnaround Times:

**Standard** (5 days in lab)

**Expedited** (1 day in lab)

- South Florida Tri-County: Cases received by 12 pm will be delivered the next day by 4 pm.
- All other Florida counties: Cases received by 12 pm will be delivered 2 days later by 12 pm.
- Nationwide: Cases received by 12 pm will be delivered 2 days later by 5 pm.

## Frequently Asked Questions

### Are these dentures repairable?

Yes. Use a normal repair bonding agent.

### Does the lab store a digital copy?

Yes! If the copy denture ever breaks, we can make a new one using the digital copy.

### Can I make changes to the teeth?

Yes, all changes can be made. This would be considered a new denture and not a copy denture. An optional try-in device can be requested.

### What tooth and pink shades are available?

Teeth: use Vita Shade guide A1-D4

Pink: Peterson TruColor shade guide (3 available colors)



**READY TO SCHEDULE? CALL OR TEXT: 561.272.6662**